


ACERBI

KIT KLX 250/300

Art. 15-94



We thank you for choosing an original ACERBIS product !

MOUNTING INSTRUCTIONS

- Remove the seat, the original gas tank, spoilers and the original pet-cock.
- Put the original pet-cock (1) on the Acerbis gas tank using the screws (2) and washers (3) provided in the mounting kit.
- Now take the original spacer (4) * from the stock tank and insert it into the Acerbis Tank (back mount).
- Place the gas tank on the bike and attach it using the washer and original screw (5)* (back mount).
- Attach the tank to the frame using the original screw (6) (side mount).
- Next attach the gas tank to the radiators using the brackets (7). You must replace the original screws (8) with the 6x25 mm screws (9) and the 6x18 mm spacers (10) provided with the kit. Attach the brackets to the bottom of the tank using the 6x12 mm bolts (11), supplied with the kit.
- Apply to the gas tank, the screw (12), provided in the kit (seat connection).
- Connect the pet-cock (1) to the carburator and remount the seat.

ATTENTION :

- * For the KLX 300 utilized the screw and washer (5) with spacer (4) provided in the kit.

ACERBIS-REAL PLASTIC FOR REAL LIFE